

MOUNTING REFERENCE GUIDE

Philips Monitoring IntelliVue Series

CIM Cable Integrated Mounting System





Content

Congred notes	3
General notes	
How to create part numbers	
Product highlights	6
Support arm pivot	
Extended pivot arm light duty	9
Extended pivot arm heavy duty	11
Height adjustable arm	13
Extended height adjustable arm	15
High position height adjustable arm	17
Extended high position height adjustable arm	19
Counter top	21
Cart light duty	23
Cart light duty kit Philips MX with FMS	25
Support arm pivot for FMS	27
Extended pivot arm light duty for FMS	29
Down post for IntelliVue FMS	31
Mount for MP2 / X2 / X3	32
Support arm pivot for MP2 / X2 / X3	33
Extended pivot arm light duty for MP2 / X2 / X3	35
Support arm pivot for MP2 / X2 / X3, without cable hook	37
Extended pivot arm light duty for MP2 / X2 / X3, without cable hook	39
Height adjustable arm for MP2 / X2 / X3, without cable hook	41
Extended height adjustable arm for MP2 / X2 / X3, without cable hook	43
Cart light duty for MP2 / X2 / X3	45
Cart light duty kit I for MP2 / X2 / X3	47
Cart light duty kit II for MP2 / X2 / X3	49
Monitor-specific mounting adapter for Philips monitoring	51
Fixings	52
Fixings for anaesthesia machines	53



General notes

Monitor Adaptations

All mounting systems for Philips IntelliVue series are prepared to accept the Philips table top mount.

Mounting solutions for the Efficia series and Sure Signs

/ Vital Signs are shown in the catalogue "5" plunge plate adaptation".

Please note!

The Philips table top mount must be ordered separately by Philips.

Components' grounding

All arms with integrated cable management are equipped with **component grounding**. No grounding is required for arms with external or semi-integrated cable management with cables not being led through joints.

Load capacity

All **height adjustable arms** have a load capacity of 22 kg / 48.5 lbs. Rear mounted monitor adaptations and horizontally

mounted 5" plunge plates accept a maximum monitor weight of 18 kg / 39.7 lbs. An additional 4 kg / 8.8 lbs may be mounted on a down post underneath the front swivel part. Depending on material and tensile factor all weight indications have a 4-fold and/ or 6-fold safety factor according to IEC 60601-1 standards.

Custom-designed solutions and MRI applications

Please contact us for **custom-designed** solutions.

Product marking

All products which are classified as medical product risk class 1 are CE marked. As **medical product risk class 1** they are conform with the regulations (EU) 2017/745 (MDR) and IEC 60601-1. Other standards of the IEC 60601 series are not applicable to our products, regardless of their version.

Warranty

All CIM arms have a 5-year warranty to be free of defects in material and fuctionality from the date of delivery.





How to create part numbers

Part numbers can be easily created by following the steps below:

Please replace the letters in green with the part no. code as per your requirements.

1. Components' grounding:

 $\mathsf{A} \qquad \qquad \mathsf{=} \qquad \mathsf{grounded}$

2. Choose the arm version required:

224PH / 230PH = pivot arm

33030PH = extended pivot arm light duty H33030PH = extended pivot arm heavy duty

431PH = height adjustable arm

73130PH = extended height adjustable arm 531PH = high position height adjustable arm

83130PH = extended high position height adjustable arm

224FMS / 230FMS = pivot arm, FMS interface

32424FMS / 33030FMS = extended pivot arm light duty, FMS interface

224MP2 / 230MP2 = pivot arm, MP2/ X2/ X3 interface

224MP2PS / 230MP2PS = pivot arm, MP2/ X2/ X3 interface with power supply enclosure

33030MP2 = extended pivot arm light duty, MP2/ X2/ X3 interface

33030MP2PS = extended pivot arm light duty, MP2/ X2/ X3 interface with power supply

224X2X3 / 230X2X3 = pivot arm, without cable hook, MP2/ X2/ X3 interface

= extended pivot arm light duty, without cable hook, MP2/ X2/ X3 interface

431X2X3 = height adjustable arm, without cable hook, MP2/ X2/ X3 interface with power supply enclosure

73130X2X3 = extended height adjustable arm, without cable hook, MP2/ X2/ X3 interface

3. Specify where you want to mount the arm and select the fixing (Y):

1V = pole \emptyset 23 - 40 mm / 0.9 - 1.6" 1 = pole \emptyset 25 mm / 1.0"

2F = horizontal rails 10 x 25, 8 x 35, 10 x 30, 10 x 50 mm

3 = vertical rail 10 x 25 mm 4S = vertical wall channel

13 = Mindray medical supply units 17 = Amico pendants & headwalls

9 Eöwenstein Medical anaestesia machine



How to create part numbers

4. Choose the colour requested for decorative parts (C):

```
gr = RAL 7024 (graphite grey) graphite = RAL 5013 (cobalt blue)
```

5. Choose the colour requested for the aluminium die-cast parts (D):

```
1 = RAL 9016 (traffic white)
2 = RAL 9002 (grey white)
```

6. Specify the spring depending on the weight of your monitor/ device - only with height adjustable arms (X):

```
0 - 3.5 \text{ kg} / 0 - 7.7 \text{ lbs (only for MP2/ X2/ X3)}
60
                                               2 - 5 kg / 4.4 - 11 lbs
                                        =
95
                                               4 - 9 kg / 8.8 - 19.8 lbs
110
                                               4.5 - 12.5 kg / 9.9 - 27.6 lbs
130
                                               7.5 - 13 kg / 16.5 - 28.7 lbs
                                        =
160
                                               12.5 - 19 kg / 27.6 - 41.9 lbs
175
                                               13 - 16 kg / 28.7 - 35.3 lbs
220
                                               16 - 22 kg / 35.3 - 48.5 lbs
                                        =
```

Example

01 Grounding		02 Arm version		03 Fixing	04 & 05 Colour		06 Spring tension
A	_	431PH	_	Υ	CD	_	X
A	_	431PH	_	4\$	gr1	_	95



Product highlights

Non-height adjustable arm

Pivot and extended pivot arms offer full cable integration and easy cleaning for infection control in all critical care areas in the clinical environment. Our 5-years warranty guarantees reliable, trouble-free and maintenance-free operation. The swivel function at the fixing point as well as the tilt/swivel function of the monitor adaptation offers optimum positioning of your device. The extended pivot arm is a two-part arm with central pivot point between the two arm sections. All bottommounted monitors have the safety feature of a 20° tilt stop.

Height adjustable arm

The height adjustable arm offers **full cable integration**, one-wipe clean surfaces and our **5-year warranty** on the entirety of the product. Vertically and horizontally adjustable, the arm offers optimal positioning of your device. A constant viewing angle when vertically adjusted is guaranteed by a parallel linkage. **Weight-optimized springs counter-balance** the monitor without making any adjustments. All bottom-mounted monitors have the safety feature of a **20° tilt stop**. Individual rotational stops may be integrated in any joint.

To provide the perfect combination of arm and monitor, the CIM articulated height adjustable arm has a simple quarter-turn lock to hold the arm safely in position when removing the monitor. Furthermore the arm is equipped with the unique safety feature of a slow-release gas spring. Should anybody try to remove the monitor before locking the arm in position, the arm will only rise at a slow and even tempo.







Support arm pivot



Monitor adaption prepared to accept Philips table top mount

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm 18 kg / 39.7 lbs 18 kg / 39.7 lbs Max. load monitor adaptation

arm type 240 **Product weight** arm type 300 1.7 kg / 3.7 lbs1.8 kg / 4.0 lbs

Part number arm type 240 arm type 300 A-224PH-Y.CD A-230PH-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp



Adapter for Löwenstein Medical anaesthesia machines







Support arm pivot

Technical data - arm type 300

SIDE VIEW




TOP VIEW







Extended pivot arm light duty



Monitor adaption prepared to accept Philips table top mount

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

 $\begin{array}{lll} \textbf{Max. load capacity support arm} & 8 \text{ kg} \text{ / } 17.6 \text{ lbs} \\ \textbf{Max. load monitor adaptation} & 8 \text{ kg} \text{ / } 17.6 \text{ lbs} \\ \textbf{Product weight} & 2.7 \text{ kg} \text{ / } 5.9 \text{ lbs} \\ \end{array}$

Part number A-33030PH-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended pivot arm light duty

Technical data

SIDE VIEW




TOP VIEW







Extended pivot arm heavy duty



Monitor adaption prepared to accept Philips table top mount

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

 $\begin{array}{lll} \textbf{Max. load capacity support arm} & 18 \ kg \ / \ 39.7 \ lbs \\ \textbf{Max. load monitor adaptation} & 18 \ kg \ / \ 39.7 \ lbs \\ \textbf{Product weight} & 3.0 \ kg \ / \ 6.6 \ lbs \\ \end{array}$

Part number A-H33030PH-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended pivot arm heavy duty

Technical data

SIDE VIEW



TOP VIEW







Height adjustable arm



Monitor adaption prepared to accept Philips table top mount

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Product weight 3.7 kg / 8.2 lbs

Balancing range 600 N (2-5 kg / 4.4-11 lbs); 950 N (4-9 kg / 8.8-19.8 lbs); 1100 N (4.5-12.5 lbs); 1100 N (4.5-12.5 lbs); 1100 lbs); 1100 lbs

9.9-27.6 lbs); 1300 N (7.5-13 kg / 16.5-28.7 lbs); 1600 N (12.5-19 kg / 27.6-41.9 lbs);

1750 N (13-16 kg / 28.7-35.3 lbs); 2200 N (16-22 kg / 35.3-48.5 lbs)

Part number A-431PH-Y.CD X

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

100 mm 10

Adapter for Löwenstein Medical anaesthesia machines







Height adjustable arm

Technical data

SIDE VIEW



TOP VIEW







Extended height adjustable arm



Monitor adaption

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm 22 kg / 48.5 lbs Max. load monitor adaptation 22 kg / 48.5 lbs **Product weight** 4.8 kg / 10.6 lbs

Balancing range 600 N (2-5 kg / 4.4-11 lbs); 950 N (4-9 kg / 8.8-19.8 lbs); 1100 N (4.5-12.5 kg /

9.9-27.6 lbs); 1300 N (7.5-13 kg / 16.5-28.7 lbs); 1600 N (12.5-19 kg / 27.6-41.9 lbs);

1750 N (13-16 kg / 28.7-35.3 lbs); 2200 N (16-22 kg / 35.3-48.5 lbs)

Part number A-73130PH-Y.CD X

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended height adjustable arm

Technical data

SIDE VIEW



TOP VIEW









Monitor adaption prepared to accept Philips table top mount

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm 22 kg / 48.5 lbs Max. load monitor adaptation 22 kg / 48.5 lbs

Product weight 4.1 kg / 9.0 lbs

Balancing range $600 \text{ N} \text{ (2-5 kg } / \text{ 4.4-11 lbs)}; 950 \text{ N} \text{ (4-9 kg } / \text{ 8.8-19.8 lbs)}; 1100 \text{ N} \text{ (4.5-12.5 kg } / \text{ 8.8-19.8 lbs)}; 1100 \text{$

9.9-27.6 lbs); 1300 N (7.5-13 kg / 16.5-28.7 lbs); 1600 N (12.5-19 kg / 27.6-41.9 lbs);

1750 N (13-16 kg / 28.7-35.3 lbs); 2200 N (16-22 kg / 35.3-48.5 lbs)

Part number A-531PH-Y.CD_X

Fixings



Adapter for vertical wall channel



Vertical rail clamp

details page 52



Universal pole clamp



1F 1C Fixed pole clamp

Adapter for Mindray medical supply units



Adapter for Amico pendants & headwalls

CIM-MED.COM

17

17







High position height adjustable arm

Technical data

SIDE VIEW



TOP VIEW









Monitor adaption prepared to accept Philips table top mount

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm 22 kg / 48.5 lbs Max. load monitor adaptation 22 kg / 48.5 lbs

Product weight 5.2 kg/ 11.5 lbs

Balancing range $600 \text{ N} \text{ (2-5 kg } / \text{ 4.4-11 lbs)}; 950 \text{ N} \text{ (4-9 kg } / \text{ 8.8-19.8 lbs)}; 1100 \text{ N} \text{ (4.5-12.5 kg } / \text{ 8.8-19.8 lbs)}; 1100 \text{$

9.9-27.6 lbs); 1300 N (7.5-13 kg / 16.5-28.7 lbs); 1600 N (12.5-19 kg / 27.6-41.9 lbs);

1750 N (13-16 kg / 28.7-35.3 lbs); 2200 N (16-22 kg / 35.3-48.5 lbs)

Part number A-83130PH-Y.CD_X

Fixings



Adapter for vertical wall channel



Vertical rail clamp

details page 52



Universal pole clamp



1C Fixed pole clamp

1F



Adapter for Mindray

medical supply units



Adapter for Amico pendants & headwalls

17

19







Extended high position height adjustable arm

Technical data

SIDE VIEW



TOP VIEW







Counter top



Monitor adaption

Colours to choose

Max. load monitor adaptation

Part number

prepared to accept Philips table top mount

decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts $\,$ RAL 9016 traffic white / RAL 9002 grey white

18 kg / 39.7 lbs

counter top with base plate (see illustration)

T-HPH-60.CD_50 height 50 mm / 2"

T-HPH-60.CD_120 height 120 mm / 4.7"

counter top without base plate

T-HPH-69.CD_50 height 50 mm / 2" T-HPH-69.CD_120 height 120 mm / 4.7"







Counter top

Technical data - with base plate

SIDE VIEW




TOP VIEW





Cart light duty



Monitor adaption prepared to accept Philips table top mount

Accessories 5-star base aluminum base with counter weight

double wheel castors, ø 100 mm/ 3.9", 2 lockable

column to accept components and accessories on opposite sides

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

base RAL 9016 traffic white

handle RAL 9016 traffic white / RAL 5013 cobalt blue / RAL 7024 graphite grey

Max. load capacity 14 kg / 30.9 lbs

Optional accessories basket big 300 x 200 x 180 mm/ 11.8 x 7.9 x 7.1"

basket small 200 x 150 x 100 mm/ 7.9 x 5.9 x 3.9"

handle

Part number cart R-002PH.CD

basket big Z-51-4Z.1 basket small Z-52-4Z.1 handle Z-HD-4.C1







Cart light duty

Technical data

SIDE VIEW




TOP VIEW









prepared to accept Philips table top mount **Monitor adaption**

Cart light duty 5-star base aluminum base with counter weight

double wheel castors, ø 100 mm / 3.9", 2 lockable

column to accept components and accessories on opposite sides

Colours decorative parts RAL 5013 cobalt blue

> RAL 9016 traffic white base handle RAL 9016 traffic white

Maximum load monitor adaptation

9.5 kg / 20.9 lbs **Maximum load FMS** 10 kg / 22 lbs **Maximum load basket** 2.5 kg / 5.5 lbs

Kit number R-002PH.bl1_286834

> Consisting of the following parts: R-002PH.bl1 cart light duty

Z-HD-4.w1 handle for carts Z-FMS-4Z.e1_56 FMS adaptation

Z-52-4Z.1 basket small 200 x 150 x 100 mm / 7.9 x 5.9 x 3.9"







Cart light duty kit Philips MX with FMS

Technical data

SIDE VIEW



TOP VIEW







Support arm pivot for FMS



Monitor adaption adapter for FMS

Colours to choosedecorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm18 kg / 39.7 lbsMax. load monitor adaptation18 kg / 39.7 lbs

Product weight arm type 240 arm type 300 1.6 kg/ 3.5 lbs 1.7 kg/ 3.8 lbs

Part number arm type 240 arm type 300
A-224FMS-Y.CD A-230FMS-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp



Adapter for Löwenstein Medical anaesthesia machines







Support arm pivot for FMS

Technical data - arm type 300

SIDE VIEW



TOP VIEW







Extended pivot arm light duty for FMS



Monitor adaption adapter for FMS

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm 8 kg / 17.6 lbs 8 kg / 17.6 lbs Max. load monitor adaptation **Product weight**

arm type 240 arm type 300 2.4 kg / 5.3 lbs2.6 kg / 5.7 lbs

Part number arm type 240 arm type 300 A-32424FMS-Y.CD A-33030FMS-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended pivot arm light duty for FMS

Technical data - arm type 300

SIDE VIEW



TOP VIEW







Down post for IntelliVue FMS





Part number

To attach to non-height adjustable arms

Z-DPR-70.CD_150 rotatable, length 150 mm / 5.9", max. load 12 kg / 26.5 lbs Z-DPR-70.CD_300 rotatable, length 300 mm / 11.8", max. load 10 kg / 22 lbs

To attach to height adjustable arms

Z-DPR-71.CD_150 rotatable, length 150 mm / 5.9", max. load 12 kg / 26.5 lbs Z-DPR-71.CD_300 rotatable, length 300 mm / 11.8", max. load 10 kg / 22 lbs





Mount for MP2/ X2/ X3

MP2/ X2/ X3 mount with power supply

Part no. Z-MP2PS-6.bl1 univeral horizonal rail clamp

Part no. Z-MP2PS-3.1 vertical rail clamp

Part no. Z-MP2PS-4Z.1 adapter for vertical wall channel (illustration)

Part no. Z-MP2PS-1V.1 universal pole clamp



MP2/ X2/ X3 mount without power supply

Part no. Z-MP2-6.bl1 univeral horizonal rail clamp

Part no. Z-MP2-3.1 vertical rail clamp

Part no. Z-MP2-4Z.1 adapter for vertical wall channel (illustration)


Part no. Z-MP2-1V.1 universal pole clamp







Support arm pivot for MP2/ X2/ X3



Monitor adaption

Colours to choose

Max. load capacity support arm Max. load monitor adaptation

Product weight

Part number

adapter for MP2/ X2/ X3

decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey aluminum parts RAL 9016 traffic white / RAL 9002 grey white

18 kg / 39.7 lbs 18 kg / 39.7 lbs

arm type 240 arm type 300 2.1 kg / 4.6 lbs2.2 kg / 4.9 lbs

without power supply (illustration) 2.3 kg / 5.1 lbs2.5 kg / 5.5 lbs with power supply

arm type 240 arm type 300 A-224MP2-Y.CD A-230MP2-Y.CD

without power supply (illustration)

A-224MP2PS-Y.CD A-230MP2PS-Y.CD with power supply

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp



Adapter for Löwenstein Medical anaesthesia machines







Support arm pivot for MP2/ X2/ X3

Technical data - arm type 300

SIDE VIEW




TOP VIEW







Extended pivot arm light duty for MP2/ X2/ X3



Monitor adaption adapter for MP2/ X2/ X3

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm8 kg / 17.6 lbsMax. load monitor adaptation8 kg / 17.6 lbs

Product weight 3.2 kg / 7.1 lbs without power supply (illustration)

3.4 kg / 7.5 lbs with power supply

Part number A-33030MP2-Y.CD without power supply (illustration)

A-33030MP2PS-Y.CD with power supply

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

9

Adapter for Löwenstein Medical anaesthesia machines







Extended pivot arm light duty for MP2/X2/X3

Technical data

SIDE VIEW



TOP VIEW







Support arm pivot for MP2/ X2/ X3, without cable hook



Monitor adaption

Colours to choose

Max. load capacity support arm Max. load monitor adaptation Product weight

Part number

adapter for MP2/ X2/ X3

decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey aluminum parts RAL 9016 traffic white / RAL 9002 grey white

18 kg / 39.7 lbs 18 kg / 39.7 lbs

 $\begin{array}{lll} \text{arm type 240} & \text{arm type 300} \\ 1.5 \text{ kg} \ / \ 3.3 \text{ lbs} & 1.7 \text{ kg} \ / \ 3.7 \text{ lbs} \\ \end{array}$

arm type 240 arm type 300 A-224X2X3-Y.CD A-230X2X3-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Support arm pivot for MP2/ X2/ X3, without cable hook

Technical data - arm type 300

SIDE VIEW




TOP VIEW







Extended pivot arm light duty for MP2/ X2/ X3, without cable hook



Monitor adaption adapter for MP2/ X2/ X3

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm8 kg / 17.6 lbsMax. load monitor adaptation8 kg / 17.6 lbs

Product weight 2.6 kg / 5.7 lbs

Part number A-33030X2X3-Y.CD

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended pivot arm light duty for MP2/ X2/ X3, without cable hook

Technical data

SIDE VIEW



TOP VIEW









Monitor adaption adapter for MP2/ X2/ X3

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Max. load capacity support arm22 kg / 48.5 lbsMax. load monitor adaptation18 kg / 39.7 lbs

Product weight 3.5 kg / 7.7 lbs

Balancing range 240 N (0-3.5 kg / 0-7.7 lbs); 600 N (2-5 kg / 4.4-11 lbs)

Part number A-431X2X3-Y.CD_X

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units



Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Height adjustable arm for MP2/ X2/ X3, without cable hook

Technical data

SIDE VIEW



TOP VIEW









Monitor adaption adapter for MP2/ X2/ X3

Height adjustment slow-flowing spring, parallel linkage for a constant viewing angle

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

aluminum parts RAL 9016 traffic white / RAL 9002 grey white

Product weight 4.5 kg / 9.9 lbs

Balancing range 240 N (0-3.5 kg / 0-7.7 lbs); 600 N (2-5 kg / 4.4-11 lbs)

Part number A-73130X2X3-Y.CD_X

Fixings



Adapter for vertical wall channel



Adapter for Mindray medical supply units

2F

Horizontal rail clamp



Adapter for Amico pendants & headwalls

details page 52 and 53



Vertical rail clamp



Universal pole clamp



Fixed pole clamp

Adapter for Löwenstein Medical anaesthesia machines







Extended height adjustable arm for MP2/ X2/ X3, without cable hook

Technical data

SIDE VIEW



TOP VIEW





Cart light duty for MP2 / X2 / X3



Monitor adaption adapter for MP2/ X2/ X3

Accessories 5-star base aluminum base with counter weight

double wheel castors, ø 100 mm / 3.9", 2 lockable

column to accept components and accessories on opposite sides

Colours to choose decorative parts RAL 5013 cobalt blue / RAL 7024 graphite grey

base RAL 9016 traffic white

handle RAL 9016 traffic white / RAL 5013 cobalt blue / RAL 7024 graphite grey

 $\textbf{Max. load capacity} \hspace{1.5cm} 7.0 \text{ kg } / \text{ 15.4 lbs}$

Optional accessories basket big 300 x 200 x 180 mm/ 11.8 x 7.9 x 7.1"

basket small 200 x 150 x 100 mm/ 7.9 x 5.9 x 3.9"

handle

Part number cart R-002X2X3.CD

basket big Z-51-4Z.1 basket small Z-52-4Z.1 handle Z-HD-4.C1







Cart light duty for MP2 / X2 / X3

Technical data

SIDE VIEW



TOP VIEW







Cart light duty kit I for MP2 / X2 / X3



Monitor adaption adapter for MP2/ X2/ X3

Cart light duty 5-star base aluminum base without counter weight

double wheel castors, ø 100 mm / 3.9", 2 lockable

column to accept components and accessories on opposite sides

Colours decorative parts RAL 5013 cobalt blue

> base RAL 9016 traffic white handle RAL 9013 cobalt blue

Maximum load monitor adaptation

2 kg / 4.4 lbs Maximum load basket 2.5 kg / 5.5 lbs

Kit number PH-002X2X3-01.bl1

> Consisting of the following parts: R-002X2X3.bl1 cart light duty

Z-HD-4.bl1 handle for carts

basket small 200 x 150 x 100 mm / 7.9 x 5.9 x 3.9" Z-52-4Z.1







Cart light duty kit I for MP2 / X2 / X3

Technical data

SIDE VIEW




TOP VIEW





Cart light duty kit II for MP2 / X2 / X3



Monitor adaption adapter for MP2/ X2/ X3

Cart light duty 5-star base aluminum base with counter weight

double wheel castors, ø 100 mm / 3.9", 2 lockable

column to accept components and accessories on opposite sides

Colours decorative parts RAL 5013 cobalt blue

> base RAL 9016 traffic white handle RAL 9013 cobalt blue

Maximum load monitor adaptation

2 kg / 4.4 lbs Maximum load basket 10 kg / 22 lbs

Kit number PH-002X2X3-02.bl1

Consisting of the following parts:

R-002X2X3.bl1 cart light duty Z-HD-4.bl1 handle for carts

Z-52-4Z.1 basket big 300 x 200 x 180 mm / 11.8 x 7.9 x 7.1"

Z-CT-7.1 counter weight







Cart light duty kit II for MP2 / X2 / X3

Technical data

SIDE VIEW



TOP VIEW





Monitor-specific mounting adapter for Philips monitoring



Monitor Description	Description	Part No.
Sure Signs / Vital Signs	Mounting adapter for 5" plunge plate	PH-AP-SS
Efficia CM100 / CM120 / CM150	Mounting adapter for 5" plunge plate	PH-AP-EF

Adapter plates are intended for CIM products only.

These part numbers are valid only when ordered in combination with CIM support arms, carts or counter tops.

Adapter plates ordered separately are subject to prior approval. Please note that in this case part numbers and prices will change.

Adapter plates for Sure Signs/ Vital Signs require a 5" plunge plate. You will find mounting solutions with 5" mounting plate in the catalogue: 5" plunge plate adaptation.



Fixings



Illustration	Description	Notes	Weight	Part no.
	Adapter for vertical wall channel	To mount to CIM wall channel (compatible with GCX channel) Please order separately: CIM wall channel, see catalogue accessoires	0.6 kg / 1.3 lbs	4\$
	Horizontal rail clamp	To mount to horizontal European rails 10 x 25, 8 x 35, 10 x 30, 10 x 50 mm	1.1 kg / 2.4 lbs	2F
	Vertical rail clamp	To mount to vertical European rail 10 x 25 mm	0.5 kg / 1.1 lbs	3
100	Universal pole clamp	V-block to mount to pole ø 23 - 40 mm / 0.9 - 1.6"	0.6 kg / 1.3 lbs	1V
	Fixed pole clamp	To mount to pole ø 25 mm / 1.0" To mount to pole ø 30 mm / 1.2" To mount to pole ø 35 mm / 1.4" To mount to pole ø 38 mm / 1.5"	0.3 kg / 0.7 lbs 0.4 kg / 0.9 lbs 0.4 kg / 0.9 lbs 0.4 kg / 0.9 lbs	1 1B 1F 1C
	Adapter medical supply units	For Mindray medical supply units	0.8 kg / 1.8 lbs	13
+ 00	Adapter pendants & headwalls	For Amico pendants & headwalls	0.6 kg / 1.3 lbs	17







Illustration	Description	Notes	Weight	Part no.		
	Adapter for Löwenstein Medical anaesthesia machines	Adapter for Leon / Leon plus	0.4 kg / 0.9 lbs	9		
f your anaesthes	sia machine does not have a s	ide channel already mounted to the machine, please ord	er the following in ad	dition to 4S:		
Adapter for Dräger	Adapter for Dräger	Dräger Primus / Julian / Cato / Fabius GS / Zeus side channel adapter L 127 mm / 5"	0.5 kg / 1.1 lbs	DR-J01		
	side channel adapter L 320 mm / 12.6"	1.0 kg / 2.2 lbs	DR-J02			
30	Adapter for GE	GE Avance / Aespire / Aestiva side channel adapter L 200 mm / 7.9"	0.5 kg / 1.1 lbs	DX-AES-0		
		Side onamine adapter £ 200 mm / 7.0	0.0 kg / 1.1 lb3	DA ALO-0		



CIM med GmbH

Euro-Industriepark Margot-Kalinke-Str. 9 80939 Munich / Germany

+49 89. 978 94 08-0

F +49 89. 978 94 08-29

M info@cim-med.com

www.cim-med.com

